Jeringa STERILE A + Aguja de seguridad STERILE EO

30°C

1°C

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CE 0318

(Hialuronato sódico 2% (20mg/ml) jeringa intraarticular)

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
DESCRIPTION: Sodium hyaluronate is a natural polymer belonging to the class of glycosaminoglycans (acid mucopolysaccharides) widely distributed in animal and human connective tissues. It is a polysaccharide made up of repeating units of the disaccharide formed by N-acetylglucosamine and sodium glucuronate linked by alternating 8-1,3 and 8-1,4 glycosidic bonds. It is a major compound of joint fluid and plays an important role in the physiological functions of that fluid. ADANT PLUS is a sterile solution with a 2% concentration of sodium hyaluronate obtained from *Streptococcus zoopelidemicus* by fermentation and subsequent purification processes. It is characterized by its viscoelasticity, hydroscopicity and lubricating properties. ADANT PLUS is available in a single-use prefilled syringe which facilitates administration, maintains sterility and prevents particle contamination during use.

INDICATIONS: ADANT PLUS is indicated to reduce pain and lack of mobility due to degenerative and traumatic changes of the knee joint and other synovial joints associated with osteoarthritis.

CONTRA-INDICATIONS: The product is contraindicated in patients with a history of hypersensitivity to the product.

PRECAUTIONS: This product may worsen local inflammation in those patients with osteoarthritis of the knee that courses with inflammation and effusion of the joint. Therefore, it is advisable to administer the product after such inflammation has been relieved.

Rarely, local pain may occur after administration of the product per se. To avoid pain after injection, it is recommended to keep the joint where the administration has been performed at rest. The product should be injected exactly into the intra-articular space, otherwise it may cause pain. This product should be administered with caution in patients with a history of hypersensitivity to other drugs and in patients with liver failure or dysfunction. The product should be administered with caution in patients with impaired physiological functions.

Use in children: Since the safety of this product has not been established in children, it should be administered with caution in these patients.

Use in pregnancy and lactation: The safety of the product has not been established in pregnant women. However, although no teratogenic effects have been detected in animals, it is recommended to be administered with caution to pregnant women or women who may be pregnant. It is recommended to discontinue breast-feeding during treatment as sodium hyaluronate has been shown to be excreted in the breast milk of treated animals.

WARNINGS: Keep out of the reach and sight of children. Do not use the product if its packaging is damaged. Discard any opened or damaged package. Administer only if the solution for injection is completely clear. Keep the prefilled syringe in the original packaging. The graduation on the syringe label should only be used as a guide. Never try to straighten bent needles but, on the contrary, discard it and replace it with a new one.

DOSAGE AND ADMINISTRATION: A single intra-articular administration is recommended. In most cases, one administration per year is sufficient. Administration should be into the intra-articular space of the affected knee joint and other affected synovial joints. Administration should be performed by a specialist with experience in intra-articular administration.

Precautions in administration: **Injection:** • Since this product is injected into the joint, the entire process should be performed under strict aseptic conditions.

• If small synovial fluid retention is observed, drain it by arthrocentesis. **Others:** • Take special care not to inject the product into the blood vessels. • Do not apply the product in the eyes. • Given the viscosity of the product, it is recommended to use 21G needles. • If the product is kept refrigerated (see storage conditions), allow it to reach room temperature before use. • Use the product immediately after opening, otherwise sterility cannot be guaranteed. The product should be used only once and always on a single patient, discarding the remaining product. Otherwise cross-infection may occur.

INCOMPATIBILITIES: Contact with disinfectants such as quaternary ammonium salts including benzalkonium chloride, or containing chlorhexidine should be avoided, as precipitation of the product may occur.

ADVERSE REACTIONS: After intra-articular administration of hyaluronic acid, local adverse effects such as pain (transient), swelling and/or effusion in the injected joint have been described. Hypersensitivity reactions, including anaphylactic reactions, anaphylactic shock and angioedema, have also been reported. In these cases, administration should be discontinued, the patient monitored and suitable treatment given. Systemic adverse effects occasionally reported after intra-articular administration of hyaluronic acid are: rash, urticaria, pruritis, erythema, pyrexia, chills, peripheral edema, headache, dizziness, nausea and muscle spasms (cramps). Any serious incident involving the product should be reported to the manufacturer Meiji Pharma Spain, S.A. and to the competent authority of the Member State in which the user is established.

INFORMATION TO BE PROVIDED TO THE PATIENT WITH AN IMPLANTED PRODUCT: Since it is considered an implant, this product includes an implant card that should be given to the patient once completed together with the instructions of the same, so that the patient has the correct identification of the product that has been implanted. Due to the characteristics of the product, it is not considered necessary to make any warning or to take any precautions or special measures in addition to those already indicated in the corresponding sections, neither by the patient nor by the health professional with respect to reciprocal interferences with external influences, medical examinations or reasonably foreseeable environmental conditions.

Once implanted, no follow-up of the product is required, beyond that required by the disease itself. The substances to which the patient is exposed once the product is implanted are described qualitatively and quantitatively in the composition section.

PRESERNTATION: ADANT PLUS is presented in a single-use glass syringe containing 4,9 ml of sodium hyaluronate (98 mg), with a Luer Lock connection system. The syringe contents are sterile, endotoxin-free and apyrogenic. Each package contains 1 syringe + 1 sterile BD needle 21G.

SELF LIFE: Do not use the product after the expiration date stated on the package.

STORAGE CONDITIONS: Keep the product between 1 and 30°C. If this temperature cannot be assured, keep it in the refrigerator. Protect from light. **FRAGILE:** Place the used syringe in a sharp disposal container immediately after use. Do not dispose of ADANT PLUS prefilled syringes in household waste.

Local regulations on the proper disposal of the sharps container should be followed.

MANUFACTURER: Meiji Pharma Spain, S.A. Avda. de Madrid, 94, 28802 Alcalá de Henares, Madrid (Spain).

DISTRIBUTED BY: See information on the carton. The Summary of Safety and Clinical Performance (SSCP) of this product is available in the European database of medical devices (Eudamed), where it is associated with its Basic UDI - DI: 8430433HA298ART96. Eudamed public website: <https://ec.europa.eu/tools/eudamed>.

Revision date: 11/2023

cases, one administration per year is sufficient. Administration should be into the intra-articular space of the affected knee joint and other affected synovial joints. Administration should be performed by a specialist with experience in intra-articular administration.

Precautions in administration: **Injection:** • Since this product is injected into the joint, the entire process should be performed under strict aseptic conditions.

• If small synovial fluid retention is observed, drain it by arthrocentesis. **Others:** • Take special care not to inject the product into the blood vessels. • Do not apply the product in the eyes. • Given the viscosity of the product, it is recommended to use 21G needles. • If the product is kept refrigerated (see storage conditions), allow it to reach room temperature before use. • Use the product immediately after opening, otherwise sterility cannot be guaranteed. The product should be used only once and always on a single patient, discarding the remaining product. Otherwise cross-infection may occur.

INCOMPATIBILITIES: Contact with disinfectants such as quaternary ammonium salts including benzalkonium chloride, or containing chlorhexidine should be avoided, as precipitation of the product may occur.

ADVERSE REACTIONS: After intra-articular administration of hyaluronic acid, local adverse effects such as pain (transient), swelling and/or effusion in the injected joint have been described. Hypersensitivity reactions, including anaphylactic reactions, anaphylactic shock and angioedema, have also been reported. In these cases, administration should be discontinued, the patient monitored and suitable treatment given. Systemic adverse effects occasionally reported after intra-articular administration of hyaluronic acid are: rash, urticaria, pruritis, erythema, pyrexia, chills, peripheral edema, headache, dizziness, nausea and muscle spasms (cramps). Any serious incident involving the product should be reported to the manufacturer Meiji Pharma Spain, S.A. and to the competent authority of the Member State in which the user is established.

INFORMATION TO BE PROVIDED TO THE PATIENT WITH AN IMPLANTED PRODUCT: Since it is considered an implant, this product includes an implant card that should be given to the patient once completed together with the instructions of the same, so that the patient has the correct identification of the product that has been implanted. Due to the characteristics of the product, it is not considered necessary to make any warning or to take any precautions or special measures in addition to those already indicated in the corresponding sections, neither by the patient nor by the health professional with respect to reciprocal interferences with external influences, medical examinations or reasonably foreseeable environmental conditions.

Once implanted, no follow-up of the product is required, beyond that required by the disease itself. The substances to which the patient is exposed once the product is implanted are described qualitatively and quantitatively in the composition section.

PRESERNTATION: ADANT PLUS is presented in a single-use glass syringe containing 4,9 ml of sodium hyaluronate (98 mg) and a Luer Lock connection system. The syringe contents are sterile, endotoxin-free and apyrogenic. Each package contains 1 syringe + 1 sterile BD needle 21G.

PERIODE DE PÉREMPTION: Do not use the product after the expiration date stated on the package.

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
DISPOSAL: Place the used syringe in a sharp disposal container immediately after use. Do not dispose of ADANT PLUS prefilled syringes in household waste.

Local regulations on the proper disposal of the sharps container should be followed.

MANUFACTURER: Meiji Pharma Spain, S.A. Avda. de Madrid, 94, 28802 Alcalá de Henares, Madrid (Spain).

DISTRIBUTED BY: See information on the carton. The Summary of Safety and Clinical Performance (SSCP) of this product is available in the European database of medical devices (Eudamed), where it is associated with its Basic UDI - DI: 8430433HA298ART96. Eudamed public website: <https://ec.europa.eu/tools/eudamed>.

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