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EN SYRINGE HANDLING INSTRUCTIONS

Before using the product, if kept in the refrigerator (see storage conditions), remove the case from the refrigerator twenty to thirty minutes before use, and follow the instructions below.

1. Remove the syringe from the blister.
2. Remove the stopper by turning it.
3. Firmly attach the BD EclipseTM needle (21G) included in the package to the syringe by squeezing and turning.
4. Pull back on the safety guard. Hold the syringe with one hand and remove the transparent cap with the other hand. Bevel upwards = safety device backwards. Activate the safety device immediately after injection. Place the thumb or index finger on the rough area of the safety device and press it forward on the needle until you hear a "click" confirming its activation. Also, visually confirm that the needle is covered. If not activated, dispose immediately in a sharps container. Activation of the protection mechanism may cause the fluid remaining in the needle after injection to splash slightly.

*Before administration, carefully disinfect the injection site.

*After administration, the syringe and needle can be discarded separately.

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Sodium chloride 8,52 mg

Disodium phosphate dodecahydrate..... 0,6 mg

Hydrochloric acid, sodium hydroxide and water for injection, q.s.

injected into the joint, the entire process should be performed under strict aseptic conditions. • Discontinue treatment if no symptoms of improvement are observed during the first 5-week cycle. • If small synovial fluid retention is observed, drain it by arthrocentesis. Others: • Take special care not to inject the product into the blood vessels. • Do not apply the product in the eyes. • Given the viscosity of the product, it is recommended to use 21-23G needles. • If the product is kept refrigerated (see storage conditions), allow it to reach room temperature before use. • Use the product immediately after opening, otherwise sterility cannot be guaranteed. The product should be used only once and always on a single patient, discarding the remaining product. Otherwise cross-infection may occur. **INCOMPATIBILITIES:** Contact with disinfectants such as quaternary ammonium salts including benzalkonium chloride, or containing chlorhexidine should be avoided, as precipitation of the product may occur. **ADVERSE REACTIONS:** Some patients may: • Rarely develop an allergic reaction, including rash, urticaria and pruritus and anaphylactic shock. In these cases, administration should be discontinued, the patient closely monitored and appropriate treatment given. • Pain and swelling at the injection site have been uncommonly described (usually transient after administration). • Patients may occasionally experience swelling at the injection site, including edema, erythema, very hot sensation and pain on palpation together with impaired joint movement. • Very rarely, the occurrence of shock has been described, in which case administration should be discontinued, the patient closely monitored and appropriate treatment provided. Any serious incident involving the product should be reported to the manufacturer Meiji Pharma Spain, S.A. and to the competent authority of the Member State in which the user is established. **INFORMATION TO BE PROVIDED TO THE PATIENT WITH AN IMPLANTED PRODUCT:** Since it is considered an implant, this product includes an implant card that should be given to the patient once completed together with the

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instructions de la même, so que le patient a la correcte identification du produit que a été implantado. Due to the characteristics of the product, it is not considered necessary to make any warning or to take any precautions or special measures in addition to those already indicated in the corresponding sections, neither by the patient nor by the health professional with respect to reciprocal interferences with external influences, medical examinations or reasonably foreseeable environmental conditions. Once implanted, no follow-up of the product is required, beyond that required by the disease itself. The substances to which the patient is exposed once the product is implanted are described qualitatively and quantitatively in the Composition section. **PRESNTATION:** ADANT is presented in single-use glass syringes containing 2,5 ml of sodium hyaluronate (25 mg), with a Luer Lock connection system. The syringe contents are sterile, endotoxin-free and pyrogenic. Each package contains 1, 3 or 5 syringes + 1, 3 or 5 sterile BD needles (21G). **SHELF LIFE:** Do not use the product after the expiration date stated on the package. **STORAGE CONDITIONS:** Keep the product between 1 and 30°C. If this temperature cannot be assured, keep in the refrigerator. Protect from light. **DISPOSAL:** Place the used syringe in a sharps disposal container immediately after use. Do not dispose of ADANT prefilled syringes in household waste. Local regulations on the proper disposal of the sharps container should be followed. **MANUFACTURER:** Meiji Pharma Spain, S.A. Avda. de Madrid, 94, 28802 Alcalá de Henares, Madrid (Spain). **DISTRIBUTED BY:** See information on the carton. The Summary of Safety and Clinical Performance (SSCP) of this product is available in the European database of medical devices (Eudamed), where it is associated with its Basic UDI - Dl: 8430433HA125ART5P. Eudamed public website: <https://ec.europa.eu/tools/eudamed>.
Revision date: 11/2023

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